



Overview

The **finesse** heralds a new era in affordable optical pump sources. Delivering 4-6 Watts at 532nm, continuous-wave, and with the diode fibre-coupled to the head; it is small, efficient and does not require water-cooling. These features, combined with noise <0.5% rms (1Hz – 5MHz) and M2 <1.1, make the finesse an ideal OEM laser or research tool.

Low Noise

The **finesse's** pump diode is housed within the power supply and fibre-coupled to the laser head. Thus, thermal effects within the head are minimised. What little heat is generated within the head is removed by conduction - there is no fan or water cooling required. Only high quality optical components are used, resulting in a noise specification of <0.5% rms.

The plot shows a 1000 hour test using a photodiode with a bandwidth of 6 MHz. During the test, the environmental temperature was varied from 18° to 30° C. The noise level was maintained below 1.5% rms during the test.

Stability

The **fpu 35** power supply is a highly intelligent and functional control unit. It allows the laser to be operated in power or current mode; in power mode the output power is stabilised to <1% using optical feedback to the laser head.

The temperature of all critical components, and of the housing itself, is regulated by PID temperature controllers, solidly maintaining all temperature-sensitive parameters within the cavity at their optimum values.

Beam Quality

The typical M2 value of the **finesse** beam is <1.1. The resulting TEM₀₀ beam has >98% fit to a Gaussian profile in both the X and Y directions, with an ellipticity of 1:<1.05.




Construction

Laser Quantum builds all lasers to a high standard, and the **finesse** is no exception.

The effects of shock impacts are minimised by the use of zero-stress mounts throughout the cavity, and the laser's feet are engineered to deform under stress, eliminating mechanical strain within the head.

Before shipment each **finesse** is subjected to rigorous quality assurance, in line with ISO9001. Every unit is N₂ purged and hermetically sealed. There follows a 100 hour burn-in under user-realistic conditions.

The **finesse** is geared towards ultrafast applications as an easily integrable pump-source. However, the high specification and quality that the femtosecond market demands make this laser suitable for all high spec research and industrial applications.



A typical 1000 hour stability test, employing APC.

Specifications

power	4, 6 W
wavelength	532 nm
spatial quality (M ²)	<1.1
beam size	2.3 mm
divergence	0.4 mrad
point stability ¹	<5 urad
power stability ²	<1.0 % rms
noise ³	<0.5 % rms
bandwidth	50 GHz
coherence length	6 mm
beam angle	<1 mrad
operating temp.	15 - 40 °C

¹Measured over 36 hours within typical temperature range 20 - 28 °C

²Test duration 8 hrs

³Test bandwidth 1 Hz - 100 MHz

Features

extremely low noise
automatic power control
fibre-fed
diffraction limited beam
zero-stress, permanently aligned cavity
hermetically sealed
diode 22,000 hrs MTBF, 2 yr warranty
RS232 control
physical shutter
large LCD display, advanced features
extended warranty available
no external, closed-loop cooling